

PRESSURE CARE

QubiCell PM 268



PM268 Full Mattress Replacement System

The PM268 Full Mattress Replacement System is a high quality pressure redistribution surface suitable for clients with medium to high risk of Pressure Injury. The combination of specially designed cells incorporated with a low pressure high flow pumping system, enables maximum pressure redistribution.

The stylish pump unit has a soft touch panel and is modern and simple to use. Features of the pump unit include 10 comfort level adjustments, incremented in weight levels for a range of client needs. Static mode with auto revert after 20 minutes, low pressure alarm indicated by an audible sound and visual LED indicator. Smart pressure detection technology starts and stops the pump as required, therefore increasing the pump life expectancy, and reducing noise. Multi therapy modes assist carers and therapists in providing appropriate pressure care and improve client compliance. The lock out feature ensures peace of mind once the unit is correctly set to individual needs. A power failure alarm and service alarm alert users or carers of the need for urgent action.

The full mattress replacement system incorporates 20 cells, with the top 3 cells providing a non alternating pillow section. The remaining cells have a combination of ventilated and non ventilated cells and have cell in cell technology, to provide additional protection for clients in case of power failure. The inner cells provide a nominal support when the pump is unable to supply air pressure due to prolonged times of power failure. The ventilated cells allow for low air loss circulation assisting in reducing perspiration and increasing client comfort. The fully enclosed zippered cover is stretchable, vapour permeable, waterproof, antimicrobial with welded seams consistent with Australian infection control guidelines.

The PM268 Full Mattress Replacement System has been designed to assist with client safety, including a clearly indicated head and foot end, CPR release valves and alarms.

The PM268 Full Mattress Replacement System assists in both prevention and healing of stage 1 & 2 pressure injuries or for clients at low to medium risk.



Zenith PM268 Full Mattress Replacement System

Mattress size 200cm (L) x 90cm (W) x 20cm (H) Inflated

Safe working load of 180kg

Registered on the Australian Register of

Therapeutic Goods No. 133718 **AX.PM268**



* All products listed are included on the Australian Register of Therapeutic Goods